

```

YLOF: BSR XCANAL2 (#n)
      BSR XLOF1
      BRA LB305

```

```

XLOF1: BSR XCANAL3 (#n)
XLOF: BSR GA74
      BGE GA76

```

sort d0.L = lof(d5)
d3.L = loc(d5)

cas mem:

```

MOVE (A1)+, d4
BSR LB95C
MOVE.L D3, D0 ← MOVE.L (A1)+, D3
RTS

```

```
GA76: BSR GA71
```

→ d0.L = loc(d5)

met d0.L^s = lof(d5) } disque
d3.L^s = loc(d5) }

```
MOVE.L D0, -(SP)
```

```
MOVEQ #2, d2
```

} d0.L = lof(d5)

```
BSR GA710
```

```
MOVE.L (SP)+, D3
```

```
MOVE.L D0, -(SP)
```

```
BSR XSEEK
```

seek d5, d3, 0 (cas fichier)

```
MOVE.L (SP)+, D0
```

```
RTS
```