


Modifier la source suivant buffer_{A3}, samf et TPINS

teste



SCMOD: BSR SACOT TBUF(a6), A3
 BSR SPT

SCMOD: BSR SAREC
 BSR XMISF

3te blancs terminan
 devant : et ;
 vérif mémoire

```

MOVE.L TDHAUT, A3
ADD     SACOT-TBUF(a6), A3
ADDQ   #1, A3
CMP.L  TDMAX, A3    ⊗
BGE    GE088        → erreur
LEA    SAMF=TBUF(a6), A0
TST.B  (A0)
RTS    XRRI
  
```

```

GE088: LEA  GE087, A0
        BRA SAERR1
  
```

"Source trop longue"